

# Task Procedure

## Startup

1. Control Room – LAV mics on
2. Server - Load LabRecorder
3. Server - Load app, put in details - configures streams
4. Server - OBS - setup
5. P1 - Load OPENVR stream, start stream
6. P1 - Load app, put in details - configures streams
7. P1 - OBS – setup
8. Listening Room – LAV mics on
9. P2 - Load OPENVR stream, start stream
10. P2 - Load app, put in details - configures streams
11. P2 - OBS - setup
12. Server - LabRecorder - check all recording boxes - RECORD
13. Server - OBS - RECORD
14. Server - Start host
15. P1 - OBS – RECORD
16. Control room – audio recorder record
17. Control room - Video cameras record
18. Control room – SYNC CLAP
19. P1 - Put player into VR to calibration wait area
20. P2 - OBS - RECORD
21. Listening room - Video cameras record
22. Listening room – SYNC CLAP
23. P2 - Put player into VR to calibration wait area
24. P1 - CONNECT
25. P2 – CONNECT

## Shutdown

1. Control room - Stop video cameras and audio recorder
2. Server - Load lab recorder STOP
3. Server - OBS - STOP
4. Server - CLOSE APP
5. P1 - STOP OPENVR steam
6. P1 - STOP audio steam
7. P1 - STOP OBS
8. Listening room - Stop video cameras
9. P2 - STOP OPENVR steam
10. P2 - STOP audio steam
11. P2 - STOP OBS